



UV TABLES

Transilluminators 8W

Model	Wavelength (nm)	Filter Size (mm)	Power (W)	Type
ECX-F26.MX V1	312	210 x 260	5 x 8	UV SuperBright
ECX-F26.M V1	312	210 x 260	5 x 8	UV
ECX-F26.C V1	254	210 x 260	5 x 8	UV
ECX-F20.L V1	365	200 x 200	5 x 8	UV
ECX-F20.M V1	312	200 x 200	5 x 8	UV
ECX-F20.C V1	254	200 x 200	5 x 8	UV
ECX-F20.BLUE V1	-	200 x 200	-	Blue Light
ECX-F20.W V1	-	No Filter	5 x 8	White Light
ECX-F15.M	312	150 x 150	4 x 8	UV
ECX-F15.C	254	150 x 150	4 x 8	UV



UV TABLES

Transilluminators 8W - Multiband

Model	Wavelength (nm)	Filter Size (mm)	Power (W)	Type
TCP-26.LMX V1	365 / 312	210 x 260	4 x 8 / 5 x 8	UV SuperBright
TCP-26.LM V1	365 / 312	210 x 260	4 x 8 / 5 x 8	UV
TCP-26.LC V1	365 / 254	210 x 260	4 x 8 / 5 x 8	UV
TCP-26.MC V1	312 / 254	210 x 260	4 x 8 / 5 x 8	UV
TCP-20.LM V1	365 / 312	200 x 200	4 x 8 / 5 x 8	UV
TCP-20.LC V1	365 / 254	200 x 200	4 x 8 / 5 x 8	UV
TCP-20.MC V1	312 / 254	200 x 200	4 x 8 / 5 x 8	UV



000 «Диаэм»

Москва
ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7(383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7(923) 303-0152
krsk@dia-m.ru

Казань
+7(843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru

